For More Secure IVF Laboratories







ELECTRONIC WITNESS SYSTEM



WHY IVFID?

What is Electronic Witness Security System?

It is an electronic Witness security system that ensures that sperm, oocyte and embryo samples that belong to patients being treated in IVF clinics are impeccably matched, and all data is verified using RFID chips in all stages of the treatment (Andrology, Opu, Denudation, ICSI, Fertizilation, Transfer, Freezing – Thawing and PGT). In addition, it saves time in the clinic, facilitates inspections, monitors real-time business processes, imposes responsibilities on personnel, provides detailed reports, triggers an alarm in case of non-compliance, provides an uninterrupted workflow and a stress-free working environment. IVFID is not just an electronic witnessing system. It also includes Lab Management Software feature. It also includes embryo tracking, Tank Management, Pregnancy result and statistics reports.

Development and Production

IVFID has been developed by the knowledge and experience within many years in IVF. It is constantly updated dynamically according to the conditions of the day and the needs of the clinics. +55 IVF clinics preferred our system.

Flexibility and Compatibility

IVFID offers a flexible structure that can be designed according to the workflow rules of the clinics. The system can work locally as well as fully integrated with the ART Management Software databases of the clinic. In addition, all clinical movements can be followed in real time with our mobile application.

Security, Authory and Audits

IVFID, can keep all transactions as encrypted, and IVFID can report all the records in detail. Authorization module enables to define user based authorization.

User Friendly Interface and Easy Operation


Follow ups and controls are easy wwith Multi-language system for laboratory staff. It can be made individual language selections based on user.

Installation

IVFID does not require any modifications to laminarflow cabinets. Drilling or etc. Installation can be done very easily and in a short time without the need for assembly processes.



How to work IVFID?



Biopsy

Cryopreservation Transfer

oocyte and sperm containers match

by using RFID chips.

Tank Control





- Easy Installation
- Mobile Check with RFID Reader
- 2D Tank mapping
- **Full Capacity**
- Central Management
- Concent Forms
- Deadline Follow-up and Notification
- -196 degree durable Label and RFID Chip

Mobile Application













DEGODE TEKNOLOJÍ HÍZMETLERÍ SAN. ve TÍC. LTD.ŞTÍ.

Sevenler Cd. No:6 EA21-1 Post Code: 34887

Sancaktepe - İstanbul - TÜRKİYE

+90 216 504 7928 / +90 532 582 8283

www.ivfid.com info@ivfid.com

Address

Phone

E-mail

f □ □ / ivfid